

For Research Use Only

JAKMIP2 Polyclonal antibody

Catalog Number: 16843-1-AP



Basic Information

Catalog Number:

16843-1-AP

Size:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10551

GenBank Accession Number:

BC017354

GeneID (NCBI):

9832

UNIPROT ID:

Q96AA8

Full Name:

janus kinase and microtubule
interacting protein 2

Calculated MW:

810 aa, 95 kDa

Observed MW:

95 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, human brain tissue, rat brain
tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 16843-1-AP (JAKMIP2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.