

For Research Use Only

# FGF2-Specific Polyclonal antibody

Catalog Number: 16872-1-AP



## Basic Information

Catalog Number:

16872-1-AP

Size:

207 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_002006

GeneID (NCBI):

2247

UNIPROT ID:

P09038

Full Name:

fibroblast growth factor 2 (basic)

Calculated MW:

31 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

Fibroblast growth factor-2 (FGF-2), also referred to as basic FGF, belongs to a family of growth factors that stimulate proliferation of cells of mesenchymal, epithelial and neuroectodermal origin. The mRNA for this gene contains multiple polyadenylation sites, and is alternatively translated from non-AUG (CUG) and AUG initiation codons, resulting in at least five different isoforms with distinct properties.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

