

For Research Use Only

FGF1-Specific Polyclonal antibody

Catalog Number: 16873-1-AP



Basic Information

Catalog Number:

16873-1-AP

Size:

170 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000800

GeneID (NCBI):

2246

UNIPROT ID:

P05230

Full Name:

fibroblast growth factor 1 (acidic)

Calculated MW:

17 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data