

For Research Use Only

# KLHDC4 Polyclonal antibody

Catalog Number: 16915-1-AP

Featured Product

1 Publications



## Basic Information

Catalog Number:

16915-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10096

GenBank Accession Number:

BC011680

GeneID (NCBI):

54758

UNIPROT ID:

Q8TBB5

Full Name:

kelch domain containing 4

Calculated MW:

520aa, 58 kDa; 489aa, 55 kDa

Observed MW:

58 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : mouse thymus tissue, human stomach tissue, Jurkat cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Masayo Kamikokura	38473243	Cancers (Basel)	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 16915-1-AP (KLHDC4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.