

For Research Use Only

# PRR7 Polyclonal antibody

Catalog Number: 16917-1-AP



## Basic Information

Catalog Number:

16917-1-AP

Size:

150 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10198

GenBank Accession Number:

BC004261

GeneID (NCBI):

80758

UNIPROT ID:

Q8TB68

Full Name:

proline rich 7 (synaptic)

Calculated MW:

31 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

## Background Information

PRR7 is found in a postsynaptic membrane complex with DLG4 and GRIN1. It is involved in modulating neural activities via interactions with the NMDA receptor and PSD-95, and PSD core formation.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

