

For Research Use Only

UBE2Z Polyclonal antibody

Catalog Number: 16928-1-AP



Basic Information

Catalog Number:

16928-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10347

GenBank Accession Number:

BC015169

GeneID (NCBI):

65264

UNIPROT ID:

Q9H832

Full Name:

ubiquitin-conjugating enzyme E2Z

Calculated MW:

28 kDa, 38 kDa

Observed MW:

45 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : human testis tissue, mouse pancreas tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 16928-1-AP (UBE2Z antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.