

For Research Use Only

MOBKL2A Polyclonal antibody

Catalog Number: 16942-1-AP

1 Publications



Basic Information

Catalog Number:

16942-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10540

GenBank Accession Number:

BC015049

GeneID (NCBI):

126308

UNIPROT ID:

Q96BX8

Full Name:

MOB1, Mps One Binder kinase activator-like 2A (yeast)

Calculated MW:

217 aa, 25 kDa

Observed MW:

25 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:3000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Positive Controls:

WB: mouse thymus tissue, Jurkat cells, mouse lymph tissue, Raji cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Zhi-Hang Huang	36401602	J Cell Mol Med	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 16942-1-AP (MOBK2A antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.