

For Research Use Only

MOBKL2A Polyclonal antibody

Catalog Number: 16942-1-AP **1 Publications**



Basic Information

Catalog Number: 16942-1-AP	GenBank Accession Number: BC015049	Purification Method: Antigen affinity purification
Size: 350 µg/ml	GeneID (NCBI): 126308	Recommended Dilutions: WB 1:1000-1:3000
Source: Rabbit	UNIPROT ID: Q96BX8	
Isotype: IgG	Full Name: MOB1, Mps One Binder kinase activator-like 2A (yeast)	
Immunogen Catalog Number: AG10540	Calculated MW: 217 aa, 25 kDa	
	Observed MW: 25 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB: mouse thymus tissue, Jurkat cells, mouse lymph tissue, Raji cells
Cited Applications: WB	
Species Specificity: human, mouse, rat	
Cited Species: mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Zhi-Hang Huang	36401602	J Cell Mol Med	WB


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 16942-1-AP (MOBK2A antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.