

For Research Use Only

N6AMT2 Polyclonal antibody

Catalog Number: 16954-1-AP



Basic Information

Catalog Number:

16954-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10589

GenBank Accession Number:

BC018091

GeneID (NCBI):

221143

UNIPROT ID:

Q8WVE0

Full Name:

N-6 adenine-specific DNA
methyltransferase 2 (putative)

Calculated MW:

214 aa, 25 kDa

Observed MW:

25 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB: Jurkat cells, K-562 cells, mouse ovary tissue,
mouse testis tissue, HEK-293 cells, mouse thymus
tissue, mouse kidney tissue, A431 cells, human liver
tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 16954-1-AP (N6AMT2 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.