

For Research Use Only

IL-2 Polyclonal antibody

Catalog Number: 17021-1-AP

3 Publications



Basic Information

Catalog Number:

17021-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10199

GenBank Accession Number:

BC070338

GeneID (NCBI):

3558

UNIPROT ID:

P60568

Full Name:

interleukin 2

Calculated MW:

153 aa, 18 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human

Cited Species:

mouse

Background Information

IL2 is a secreted cytokine involved in the operation and regulation of immune system. IL2 plays a pivotal role in T-cell activation and effector functions, including T cell proliferation, interferon- γ production and cytotoxicity. IL2 can stimulate the propagation of lymphocytes and induce the cytotoxic T lymphocytes and lymphokine-activated killer cells for tumor cells. IL2 also influences memory T cell homeostasis through the regulation of memory T cell numbers and enhances major histocompatibility complex II (MHC II) expression on tumor cells.

Notable Publications

Author	Pubmed ID	Journal	Application
Zhongci Hang	35180124	J Alzheimers Dis	WB
Ziqiang Wang	36683735	Regen Biomater	IF
Zhu Wang	37253929	Int Urol Nephrol	IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

