

# SOHLH2 Polyclonal antibody

Catalog Number: 17057-1-AP

## Basic Information

**Catalog Number:**

17057-1-AP

**Size:**

140 µg/ml

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG10768

**GenBank Accession Number:**

BC025383

**GeneID (NCBI):**

54937

**UNIPROT ID:**

Q9NX45

**Full Name:**

spermatogenesis and oogenesis  
specific basic helix-loop-helix 2

**Calculated MW:**

225 aa, 25 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

ELISA

**Species Specificity:**

human

## Background Information

SOHLH2, also named as TEB1, is a germ cell-specific bHLH transcription factor which may be involved in spermatogenesis and oogenesis. SOHLH2 has three isoforms with MW 25; 47 and 56 kDa.

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

