

For Research Use Only

PMFBP1 Polyclonal antibody

Catalog Number: 17061-1-AP **4 Publications**



Basic Information

Catalog Number:

17061-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10789

GenBank Accession Number:

BC028423

GeneID (NCBI):

83449

UNIPROT ID:

Q8TBY8

Full Name:

polyamine modulated factor 1
binding protein 1

Calculated MW:

1022 aa, 119 kDa

Observed MW:

119 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse

Cited Species:

human, mouse

Positive Controls:

WB : mouse testis tissue, HepG2 cells, MCF-7 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yan-Wei Sha	30298696	Clin Genet	IF
Leslie A Sandusky-Beltran	31253191	Alzheimers Res Ther	WB
Hua Nie	35860846	Mol Genet Genomic Med	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 17061-1-AP (PMFBP1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.