

For Research Use Only

DNAJB8 Polyclonal antibody

Catalog Number: 17071-1-AP **2 Publications**



Basic Information

Catalog Number:

17071-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10899

GenBank Accession Number:

BC029521

GeneID (NCBI):

165721

UNIPROT ID:

Q8NHS0

Full Name:

DnaJ (Hsp40) homolog, subfamily B, member 8

Calculated MW:

232 aa, 26 kDa

Observed MW:

26 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Positive Controls:

WB: HeLa cells, human brain tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Qianxing Zou	34825148	iScience	WB, IF
Qianxing Zou	35600921	STAR Protoc	WB, IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 17071-1-AP (DNAJB8 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.