

For Research Use Only

RAI1 Polyclonal antibody

Catalog Number: 17084-1-AP

1 Publications



Basic Information

Catalog Number:

17084-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9046

GenBank Accession Number:

BC021209

GeneID (NCBI):

10743

UNIPROT ID:

Q7Z5J4

Full Name:

retinoic acid induced 1

Calculated MW:

1906aa, 203 kDa; 966aa, 104 kDa

Observed MW:

250-260 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : HeLa cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Rong Yu	36334618	Gene	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 17084-1-AP (RAI1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.