

For Research Use Only

RUNDC3B Polyclonal antibody, PBS Only

Catalog Number: 17113-1-PBS



Basic Information

Catalog Number:

17113-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10797

GenBank Accession Number:

BC022520

GeneID (NCBI):

154661

UNIPROT ID:

Q96NLO

Full Name:

RUN domain containing 3B

Calculated MW:

387 aa, 43 kDa

Observed MW:

45 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 17113-1-AP (RUNDC3B antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17113-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed MCF-7 cells using 17113-1-AP (RUNDC3B antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 17113-1-PBS in a different storage buffer formulation.