

For Research Use Only

# GMCL1L Polyclonal antibody

Catalog Number: 17129-1-AP



## Basic Information

Catalog Number:

17129-1-AP

Size:

110 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10873

GenBank Accession Number:

BC024185

GeneID (NCBI):

64396

Full Name:

germ cell-less homolog 1  
(Drosophila)-like

Calculated MW:

526 aa, 60 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

## Background Information

GMCL1L, also named as GCL and GMCL1P1, is possible function in spermatogenesis. It enhances the degradation of MDM2 and increases the amount of p53 probably by modulating the nucleocytoplasmic transport.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data