

For Research Use Only

# ACTRT1 Polyclonal antibody, PBS Only

Catalog Number:17130-1-PBS



## Basic Information

Catalog Number:

17130-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10874

GenBank Accession Number:

BC014597

GeneID (NCBI):

139741

UNIPROT ID:

Q8TDG2

Full Name:

actin-related protein T1

Calculated MW:

376 aa, 42 kDa

Observed MW:

40-45 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

IP, Indirect ELISA

Species Specificity:

human, mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



IP result of anti-ACTRT1 (IP:17130-1-AP, 4ug; Detection:17130-1-AP 1:800) with mouse testis tissue lysate 1800 ug. This data was developed using the same antibody clone with 17130-1-PBS in a different storage buffer formulation.