

For Research Use Only

GALNTL5 Polyclonal antibody, PBS Only

Catalog Number: 17131-1-PBS



Basic Information

Catalog Number:

17131-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10876

GenBank Accession Number:

BC022021

GeneID (NCBI):

168391

UNIPROT ID:

Q7Z4T8

Full Name:

UDP-N-acetyl-alpha-D-galactosamine:polypeptide N-acetylgalactosaminyltransferase-like 5

Calculated MW:

443 aa, 51 kDa

Observed MW:

51 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 17131-1-AP (GALNTL5 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17131-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human lung using 17131-1-AP (GALNTL5 antibody) at dilution of 1:100 (under 40x lens). This data was developed using the same antibody clone with 17131-1-PBS in a different storage buffer formulation.