

For Research Use Only

MRPS33 Polyclonal antibody

Catalog Number: 17162-1-AP



Basic Information

Catalog Number:

17162-1-AP

Concentration:

500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10566

GenBank Accession Number:

BC015462

GeneID (NCBI):

51650

UNIPROT ID:

Q9Y291

Full Name:

mitochondrial ribosomal protein S33

Calculated MW:

106 aa, 13 kDa

Observed MW:

13-15 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17162-1-AP (MRPS33 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.