

For Research Use Only

# FUT11 Polyclonal antibody, PBS Only

Catalog Number:17175-1-PBS



## Basic Information

Catalog Number:

17175-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10815

GenBank Accession Number:

BC036037

GeneID (NCBI):

170384

UNIPROT ID:

Q495W5

Full Name:

fucosyltransferase 11 (alpha (1,3)  
fucosyltransferase)

Calculated MW:

492 aa, 56 kDa

Observed MW:

56 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17175-1-AP (FUT11 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17175-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human testis using 17175-1-AP (FUT11 antibody) at dilution of 1:100 (under 10x lens). This data was developed using the same antibody clone with 17175-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human testis using 17175-1-AP (FUT11 antibody) at dilution of 1:100 (under 40x lens). This data was developed using the same antibody clone with 17175-1-PBS in a different storage buffer formulation.