

For Research Use Only

LETM2 Polyclonal antibody, PBS Only

Catalog Number: 17180-1-PBS

Featured Product



Basic Information

Catalog Number:

17180-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10900

GenBank Accession Number:

BC029541

GeneID (NCBI):

137994

UNIPROT ID:

Q2VYF4

Full Name:

leucine zipper-EF-hand containing transmembrane protein 2

Calculated MW:

277 aa, 31 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using 17180-1-AP (LETM2 Antibody) at dilution of 1:50 (under 10x lens). This data was developed using the same antibody clone with 17180-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using 17180-1-AP (LETM2 Antibody) at dilution of 1:50 (under 40x lens). This data was developed using the same antibody clone with 17180-1-PBS in a different storage buffer formulation.