

For Research Use Only

WDFY2 Polyclonal antibody, PBS Only

Catalog Number:17197-1-PBS



Basic Information

Catalog Number:

17197-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10415

GenBank Accession Number:

BC014004

GeneID (NCBI):

115825

UNIPROT ID:

Q96P53

Full Name:

WD repeat and FYVE domain containing 2

Calculated MW:

400 aa, 45 kDa

Observed MW:

45 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse

Background Information

WDFY2, also named Prof, is a protein that contains WD repeats and an FYVE domain. WDFY2 regulates the membrane-bound matrix metalloproteinase MT1-MMP (MMP14) efflux by controlling endosome sorting. WDFY2 may act as a tumor suppressor by limiting VAMP3-dependent recycling to inhibit matrix metalloproteinase secretion and cell invasion (PMID: 31253801). WDFY2 in the cell nucleus directly participates in homologous recombination-mediated DNA repair (PMID: 41196680).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17197-1-AP (WDFY2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17197-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using WDFY2 antibody (17197-1-AP) at dilution of 1:200 and CoraLite®594-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-4). This data was developed using the same antibody clone with 17197-1-PBS in a different storage buffer formulation.