

For Research Use Only

MRPL17 Polyclonal antibody

Catalog Number: 17214-1-AP

1 Publications



Basic Information

Catalog Number:

17214-1-AP

Concentration:

210 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10997

GenBank Accession Number:

BC012306

GeneID (NCBI):

63875

UNIPROT ID:

Q9NRX2

Full Name:

mitochondrial ribosomal protein L17

Calculated MW:

175 aa, 20 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Elena Lavdovskaia	38992089	Nat Struct Mol Biol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

