

For Research Use Only

# ARL17 Polyclonal antibody

Catalog Number: 17231-1-AP



## Basic Information

Catalog Number: 17231-1-AP	GenBank Accession Number: BC030570	Purification Method: Antigen affinity purification
Size: N/A	GeneID (NCBI): 641522	
Source: Rabbit	UNIPROT ID: Q8IVW1	
Isotype: IgG	Full Name: ADP-ribosylation factor-like 17	
Immunogen Catalog Number: AG11072	Calculated MW: 160 aa, 17 kDa	

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data