

For Research Use Only

DUS1L Polyclonal antibody

Catalog Number: 17253-1-AP

1 Publications



Basic Information

Catalog Number:

17253-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10668

GenBank Accession Number:

BC017081

GeneID (NCBI):

64118

UNIPROT ID:

Q6P1R4

Full Name:

dihydrouridine synthase 1-like (S. cerevisiae)

Calculated MW:

473 aa, 53 kDa

Observed MW:

53 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : Jurkat cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Nathan J Yu	37965204	bioRxiv	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 17253-1-AP (DUS1L antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.