

For Research Use Only

NDUFA3 Polyclonal antibody, PBS Only

Catalog Number:17257-1-PBS



Basic Information

Catalog Number:

17257-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10912

GenBank Accession Number:

BC011021

GeneID (NCBI):

4696

UNIPROT ID:

O95167

Full Name:

NADH dehydrogenase (ubiquinone) 1
alpha subcomplex, 3, 9kDa

Calculated MW:

65aa,7 kDa; 84aa,9 kDa

Observed MW:

9 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human brain tissue were subjected to SDS PAGE followed by western blot with 17257-1-AP (NDUFA3 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17257-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 17257-1-AP (NDUFA3 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Sodium Citrate buffer (pH 6.0). This data was developed using the same antibody clone with 17257-1-PBS in a different storage buffer formulation.