

For Research Use Only

PUS7L Polyclonal antibody

Catalog Number: 17273-1-AP



Basic Information

Catalog Number:

17273-1-AP

Concentration:

450 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11030

GenBank Accession Number:

BC033621

GeneID (NCBI):

83448

UNIPROT ID:

Q9H0K6

Full Name:

pseudouridylate synthase 7 homolog
(S. cerevisiae)-like

Calculated MW:

701 aa, 81 kDa

Observed MW:

81, 45 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : DU 145 cells, HEK-293 cells, HeLa cells, PC-3 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17273-1-AP (PUS7L antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.