

UBE2DNL Polyclonal antibody

Catalog Number: 17287-1-AP

Basic Information

Catalog Number:

17287-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11132

GenBank Accession Number:

BC040290

GeneID (NCBI):

100131816

Full Name:

ubiquitin-conjugating enzyme E2D N-terminal like pseudogene

Calculated MW:

210 aa, 23 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:20-1:200

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Note-IHC: suggested antigen retrieval with **TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

Positive Controls:

IHC : human testis tissue, human brain tissue, human kidney tissue, human lung tissue, human ovary tissue, human placenta tissue, human skin tissue, human spleen tissue

Background Information

Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:


PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human testis using 17287-1-AP (UBE2DNL antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human testis using 17287-1-AP (UBE2DNL antibody) at dilution of 1:50 (under 40x lens).