

For Research Use Only

TRMT10A Polyclonal antibody

Catalog Number: 17294-1-AP



Basic Information

Catalog Number:

17294-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11172

GenBank Accession Number:

BC028373

GeneID (NCBI):

93587

UNIPROT ID:

Q8TBZ6

Full Name:

RNA (guanine-9-) methyltransferase domain containing 2

Calculated MW:

339 aa, 40 kDa

Observed MW:

40 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, HepG2 cells, Jurkat cells, NIH/3T3 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17294-1-AP (TRMT10A antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.