

For Research Use Only

HLA-DQB2 Polyclonal antibody

Catalog Number: 17315-1-AP



Basic Information

Catalog Number:

17315-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11284

GenBank Accession Number:

BC031995

GeneID (NCBI):

3120

Full Name:

major histocompatibility complex,
class II, DQ beta 2

Calculated MW:

231aa, 27 kDa; 275aa, 31 kDa

Observed MW:

27 kDa, 31 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : Raji cells, HepG2 cells, K-562 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Raji cells were subjected to SDS PAGE followed by western blot with 17315-1-AP (HLA-DQB2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.