

For Research Use Only

# HLA-DQB2 Polyclonal antibody

Catalog Number: 17315-1-AP



## Basic Information

Catalog Number:

17315-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11284

GenBank Accession Number:

BC031995

GeneID (NCBI):

3120

Full Name:

major histocompatibility complex,  
class II, DQ beta 2

Calculated MW:

231aa, 27 kDa; 275aa, 31 kDa

Observed MW:

27 kDa, 31 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB: Raji cells, HepG2 cells, K-562 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Raji cells were subjected to SDS PAGE followed by western blot with 17315-1-AP (HLA-DQB2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.