

For Research Use Only

# TSGA13 Polyclonal antibody

Catalog Number: 17326-1-AP



## Basic Information

Catalog Number: 17326-1-AP	GenBank Accession Number: BC038415	Purification Method: Antigen affinity purification
Size: 173 µg/ml	GeneID (NCBI): 114960	Recommended Dilutions: IHC 1:20-1:200
Source: Rabbit	UNIPROT ID: Q96PP4	
Isotype: IgG	Full Name: testis specific, 13	
Immunogen Catalog Number: AG11305	Calculated MW: 275 aa, 32 kDa	

## Applications

Tested Applications: IHC, ELISA	Positive Controls: IHC : human testis tissue,
Species Specificity: human	

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 17326-1-AP (TSGA13 Antibody) at dilution of 1:50 (under 10x lens).