

For Research Use Only

# SH3YL1 Polyclonal antibody

Catalog Number: 17359-1-AP



## Basic Information

Catalog Number:

17359-1-AP

Size:

270 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10641

GenBank Accession Number:

BC030778

GeneID (NCBI):

26751

UNIPROT ID:

Q96HL8

Full Name:

SH3 domain containing, Ysc84-like 1  
(S. cerevisiae)

Calculated MW:

114 aa, 12, 342 aa, 37 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

