

For Research Use Only

# LY6D Polyclonal antibody, PBS Only

Catalog Number: 17361-1-PBS

Featured Product



## Basic Information

Catalog Number:

17361-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10714

GenBank Accession Number:

BC031330

GeneID (NCBI):

8581

UNIPROT ID:

Q14210

Full Name:

lymphocyte antigen 6 complex, locus D

Calculated MW:

128 aa, 13 kDa

Observed MW:

13-14 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

## Background Information

LY6D, also named as E48 and ThB, may act as a specification marker at earliest stage specification of lymphocytes between B- and T-cell development. It marks the earliest stage of B-cell specification.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data




A431 cells were subjected to SDS PAGE followed by western blot with 17361-1-AP (LY6D antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17361-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 17361-1-AP (LY6D antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 17361-1-PBS in a different storage buffer formulation.



WB result of LY6D antibody (17361-1-AP; 1:1500; incubated at room temperature for 1.5 hours) with sh-Control and sh-LY6D transfected A431 cells. This data was developed using the same antibody clone with 17361-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human skin cancer tissue slide using 17361-1-AP (LY6D antibody) at dilution of 1:200 (under 10x lens). This data was developed using the same antibody clone with 17361-1-PBS in a different storage buffer formulation.