

For Research Use Only

# TADA2B Polyclonal antibody, PBS Only

Catalog Number: 17367-1-PBS

Featured Product



## Basic Information

Catalog Number:

17367-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10783

GenBank Accession Number:

BC013164

GeneID (NCBI):

93624

UNIPROT ID:

Q86TJ2

Full Name:

transcriptional adaptor 2 (ADA2 homolog, yeast)-beta

Calculated MW:

420 aa, 48 kDa

Observed MW:

55-62 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17367-1-AP (TADA2B antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17367-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human ovary cancer tissue slide using 17367-1-AP (TADA2B antibody) at dilution of 1:400 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 17367-1-PBS in a different storage buffer formulation.