

For Research Use Only

# GLIPR1L2 Polyclonal antibody

Catalog Number: 17369-1-AP



## Basic Information

|                                      |   |   |
|--------------------------------------|---|---|
| Catalog Number:<br>17369-1-AP        | GenBank Accession Number:<br>BC029557           | Purification Method:<br>Antigen affinity purification |
| Size:<br>N/A                         | GeneID (NCBI):<br>144321                        |   |
| Source:<br>Rabbit                    | UNIPROT ID:<br>Q4G1C9                           |   |
| Isotype:<br>IgG                      | Full Name:<br>GLI pathogenesis-related 1 like 2 |   |
| Immunogen Catalog Number:<br>AG10859 | Calculated MW:<br>253 aa, 29 kDa                |   |

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

