

For Research Use Only

# GPBP1L1 Polyclonal antibody

Catalog Number: 17375-1-AP



## Basic Information

Catalog Number:

17375-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11180

GenBank Accession Number:

BC022218

GeneID (NCBI):

60313

UNIPROT ID:

Q9HC44

Full Name:

GC-rich promoter binding protein 1-like 1

Calculated MW:

474 aa, 52 kDa

Observed MW:

58 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:5000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, Jurkat cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17375-1-AP (GPBP1L1 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.