

For Research Use Only

PRPF40A Polyclonal antibody

Catalog Number: 17392-1-AP



Basic Information

Catalog Number:

17392-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11290

GenBank Accession Number:

BC027178

GeneID (NCBI):

55660

UNIPROT ID:

O75400

Full Name:

PRP40 pre-mRNA processing factor 40 homolog A (S. cerevisiae)

Calculated MW:

957 aa, 109 kDa

Observed MW:

130-140 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HEK-293 cells, HeLa cells, mouse brain tissue

IP : HEK-293 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 17392-1-AP (PRPF40A Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.

IP result of anti-PRPF40A (IP:17392-1-AP, 4ug; Detection:17392-1-AP 1:500) with HEK-293 cells lysate 3000ug.