

For Research Use Only

SLC9A11 Polyclonal antibody

Catalog Number: 17398-1-AP



Basic Information

Catalog Number:

17398-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11351

GenBank Accession Number:

BC042592

GeneID (NCBI):

284525

UNIPROT ID:

Q5TAH2

Full Name:

solute carrier family 9, member 11

Calculated MW:

1124 aa, 129 kDa

Observed MW:

135-140 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse liver tissue, HepG2 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 17398-1-AP (SLC9A11 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.