

For Research Use Only

TESK2 Polyclonal antibody

Catalog Number: 17414-1-AP



Basic Information

Catalog Number:

17414-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10964

GenBank Accession Number:

BC033085

GeneID (NCBI):

10420

UNIPROT ID:

Q96S53

Full Name:

testis-specific kinase 2

Calculated MW:

542 aa, 60 kDa

Observed MW:

62 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : NCCIT cell,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



NCCIT cell were subjected to SDS PAGE followed by western blot with 17414-1-AP (TESK2 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.