

For Research Use Only

Noxa Polyclonal antibody

Catalog Number: 17418-1-AP

2 Publications



Basic Information

Catalog Number:

17418-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11387

GenBank Accession Number:

BC032663

GeneID (NCBI):

5366

UNIPROT ID:

Q13794

Full Name:

phorbol-12-myristate-13-acetate-induced protein 1

Calculated MW:

136 aa, 15 kDa

Observed MW:

6-15 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : 5637 cells, HEK-293 cells, Karpas-299 cells, Raji cells, SCaBER cells, U-251 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Fangbing Liu	31115744	Target Oncol	WB
Huiying Zhao	28573788	Cancer Med	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17418-1-AP (Noxa antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.