

For Research Use Only

Noxa Polyclonal antibody, PBS Only

Catalog Number:17418-1-PBS



Basic Information

Catalog Number:

17418-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11387

GenBank Accession Number:

BC032663

GeneID (NCBI):

5366

UNIPROT ID:

Q13794

Full Name:

phorbol-12-myristate-13-acetate-induced protein 1

Calculated MW:

136 aa, 15 kDa

Observed MW:

6-15 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17418-1-AP (Noxa antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 17418-1-PBS in a different storage buffer formulation.