

For Research Use Only

CSGALNACT2 Polyclonal antibody

Catalog Number: 17420-1-AP **1 Publications**



Basic Information

Catalog Number:

17420-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11436

GenBank Accession Number:

BC030268

GeneID (NCBI):

55454

UNIPROT ID:

Q8N6G5

Full Name:

chondroitin sulfate N-acetylgalactosaminyltransferase 2

Calculated MW:

542 aa, 63 kDa

Observed MW:

63-70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Positive Controls:

WB : K-562 cells, RAW264.7

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Margery L Evans	29886539	Methods Mol Biol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cells were subjected to SDS PAGE followed by western blot with 17420-1-AP (CSGALNACT2 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.