

For Research Use Only

# GALNT11 Polyclonal antibody

Catalog Number: 17429-1-AP



## Basic Information

Catalog Number:

17429-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11495

GenBank Accession Number:

BC059377

GeneID (NCBI):

63917

UNIPROT ID:

Q8NCW6

Full Name:

UDP-N-acetyl-alpha-D-galactosamine:polypeptide N-acetyl(galactosaminyltransferase 11 (GalNAc-T11)

Calculated MW:

608 aa, 69 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

