

For Research Use Only

LYZL2 Polyclonal antibody

Catalog Number: 17433-1-AP **1 Publications**



Basic Information

Catalog Number:

17433-1-AP

Concentration:

230 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11502

GenBank Accession Number:

BC066294

GeneID (NCBI):

119180

UNIPROT ID:

Q7Z4W2

Full Name:

lysozyme-like 2

Calculated MW:

194 aa, 22 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

