

For Research Use Only

LYZL4 Polyclonal antibody

Catalog Number: 17443-1-AP **1 Publications**



Basic Information

Catalog Number: 17443-1-AP	GenBank Accession Number: BC016747	Purification Method: Antigen affinity purification
Size: 700 µg/ml	GeneID (NCBI): 131375	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	UNIPROT ID: Q96KX0	
Isotype: IgG	Full Name: lysozyme-like 4	
Immunogen Catalog Number: AG11532	Calculated MW: 146 aa, 16 kDa	
	Observed MW: 16 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse testis tissue, Raji cells
Cited Applications: WB, CoIP	
Species Specificity: human, mouse, rat	
Cited Species: mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Sara Bilekova	37691832	Front Cell Dev Biol	CoIP, WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 17443-1-AP (LYZL4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.