

For Research Use Only

LYZL4 Polyclonal antibody

Catalog Number: 17443-1-AP

1 Publications



Basic Information

Catalog Number:

17443-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11532

GenBank Accession Number:

BC016747

GeneID (NCBI):

131375

UNIPROT ID:

Q96KX0

Full Name:

lysozyme-like 4

Calculated MW:

146 aa, 16 kDa

Observed MW:

16 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, CoIP

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Positive Controls:

WB : mouse testis tissue, Raji cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Sara Bilekova	37691832	Front Cell Dev Biol	CoIP, WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 17443-1-AP (LYZL4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.