

For Research Use Only

TBC1D21 Polyclonal antibody

Catalog Number: 17463-1-AP



Basic Information

Catalog Number: 17463-1-AP	GenBank Accession Number: BC033516	Purification Method: Antigen affinity purification
Size: 500 µg/ml	GeneID (NCBI): 161514	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q8IYX1	
Isotype: IgG	Full Name: TBC1 domain family, member 21	
Immunogen Catalog Number: AG11301	Calculated MW: 336 aa, 39 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse testis tissue, rat testis tissue
Species Specificity: human, mouse, rat	

Background Information

TBC1D21 genes (also known as male germ cell RAB GTPase-activating protein MGCRAAGAP) are related to sterility, as determined through cDNA microarray testing of human testicular tissues exhibiting spermatogenic defects. TBC1D21 is a protein specifically expressed in the testes that exhibits specific localizations of elongating and elongated spermatids during mammalian spermiogenesis.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17463-1-AP (TBC1D21 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.