

For Research Use Only

# BDNF Polyclonal antibody

Catalog Number: 17465-1-AP **7 Publications**



## Basic Information

<b>Catalog Number:</b> 17465-1-AP	<b>GenBank Accession Number:</b> BC029795	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 310 µg/ml	<b>GeneID (NCBI):</b> 627	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P23560	
<b>Isotype:</b> IgG	<b>Full Name:</b> brain-derived neurotrophic factor	
<b>Immunogen Catalog Number:</b> AG11329	<b>Calculated MW:</b> 247 aa, 28 kDa	

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
WB

**Species Specificity:**  
human, mouse, rat

**Cited Species:**  
rat, mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Zhiyuan Wang	36438503	Theranostics	WB
Qiaozhen Qin	35883153	Stem Cell Res Ther	WB
Sen Yin	32802036	Neural Plast	WB

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data