

For Research Use Only

# BDNF Polyclonal antibody

Catalog Number: 17465-1-AP **13 Publications**



## Basic Information

Catalog Number:

17465-1-AP

Concentration:

307 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11329

GenBank Accession Number:

BC029795

GeneID (NCBI):

627

UNIPROT ID:

P23560

Full Name:

brain-derived neurotrophic factor

Calculated MW:

247 aa, 28 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

mouse, rat

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Zhiyuan Wang	36438503	Theranostics	WB
Qiaozhen Qin	35883153	Stem Cell Res Ther	WB
Sen Yin	32802036	Neural Plast	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

